健康快拍 Health Gather

Give you 24 hours of silent companionship

Principles of Pathology

In anemia diagnosis, the color characteristics of the eyelid, particularly the degree of pallor, are key indicators. In healthy individuals, the eyelid appears reddish, reflecting sufficient levels of hemoglobin in the blood. However, when anemia occurs, hemoglobin levels drop, and the color of the eyelid tends to become pale, as the blood's oxygen-carrying capacity decreases, affecting oxygen supply to tissues.

Technically, changes in eyelid color can be captured using image processing methods. Specifically, the image processing system can decompose the eyelid image into red, green, and blue channels, and compare the value of the red channel against other color channels. Studies have shown that anemic patients often exhibit lower red channel values, with RGB channel values being similar, causing the eyelid to appear pale or pinkish. This phenomenon serves as a potential indicator for quantifying hemoglobin levels, providing a non-invasive approach for evaluating anemia.

Detection Process

Selecting a picture from your gallery or taking a new photo with the device's camera

Done

Use our app to crop

the eye area

Image-based symptom detection

Predict result

Predict Hb value is 11.48

Application Principles

Data Preprocessing

Deep Learning Classification

To address the issue of glare in the images, we implemented two methods. First, we used the HSV color space to filter out glare points with excessively high brightness or low saturation. The second method involved converting the image to grayscale and using threshold operations to detect and correct over-bright areas.

During the image segmentation stage, the UNet model was employed to automatically segment the eyelid region, allowing for precise extraction of anemia-related features. The UNet architecture consists of an encoder and a decoder: the encoder progressively captures critical features within the image to locate the eyelid, while the decoder restores these features to accurately mark the eyelid area. In deep learning classification, a deep learning model is employed to classify anemia images, using a hemoglobin (Hb) concentration threshold of 12 g/dL to determine anemia presence. The model leverages multi-level feature extraction to automatically identify anemia-related features, such as color variations in the eyelid area. To improve the model's generalization, data augmentation techniques are applied during training to diversify the dataset, enhancing the model's performance on new data and strengthening the stability and accuracy of the classification results. The model achieved an accuracy of 85% to 90% on the test set, demonstrating its reliability in anemia detection.

Regression Analysis

Regression analysis is used to quantify the severity of a patient's anemia by examining trends in hemoglobin levels. This analysis selects features closely related to eyelid color—such as hue, saturation, brightness, the R channel, the difference between R and GB channels, and grayscale—to estimate hemoglobin levels. Linear and polynomial regression models are employed to precisely predict hemoglobin fluctuations, while mean squared error (MSE) and the coefficient of determination (R²) are used to assess model performance, ensuring both predictive accuracy and interpretability.

Clinical Trials

Our team, having received approval from the Institutional Review Board, has initiated clinical trials at Shuang Ho Hospital and the Hemodialysis Center of New Civil Hospital.

<u>Taipei Medical University</u> Certificate of TMU-JIRB Approval

TAULHIB No. NNR201117 Princed Tene Step et a pinnet junctic and scenia with a marphone Principal research principal and pinnet and pinnet pinnet. Solid Step 2011; Step 201; Step 2011; Step 2011; Step 2011; Step 2011; Step 2011;

The above study will be approve by expedited review process of the TMU-Joint Institutional Review Board in meeting #111-07-2(date:2022/07/12), duration of validity is from 2022/06/18 to 2023/06/17, and must be monitored by TMU-JIRB.

According to Ministry of Health and Welfare and the relevant regulations, follow-up procedures and requirements are as below: 1. Continuing Report: Continuous report frequency is every 12 months. The report shuould be

 Community report: Communities report impairing its very is mounter the report innovation of submitted in 2 months before the end of validity (2023-04-17). The trial/study cannot going if the continuous report not approve yet.
 Final Report: The report isloadd be submitted when the trial/study complete. TMU-JIRB will

2. ritilia appet: rate ruport data statisticative for supervision stores and statisticative for an experimental statisticative for supervision stores and statisticative for supervision stores and statisticative for supervision statisticative s

Chairman

4) the W. Shen 42

TMU-IIRR Form076/2020031

laune Date: 2022/06/18

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+ Φ104,0 R,1 9 20 0 4/R d.R.
The TBE-Joint Institutional Review Board performs its functions according to written
 operating procedures and couplies with 602 and with the applicable regulatory requirements.
 199bc05056a16dac22030c8
 TMU-BRB Feem076/2020031

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TMI-JIRE Form(64/20200317

Phase Results

- To date, we have successfully enrolled over 350 participants, and recruitment is actively ongoing.
- A preliminary design for the measurement mobile application has been completed.
- An anemia detection model was initially developed with an accuracy of over 94%.

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2024/06 InnoVex

Call for support our projects

We have secured for multiple patents for the relevant technologies to further ensure the legality and commercial applicability of our research findings, while protecting our intellectual property rights under the law.

uspto 🕅	TED STATES Ent and trademark off	Page 1 P.O. Box Alexandria, VA 22313 www.uspb				
ELE		NOWL		RECEIPT		
APPLICATION # 18/398,320	RECEIPT DATE / TIME 04/16/2024 08:04:37 AM Z ET			ATTORNEY DOCKET # 68507-330		
Title of Invention	0					
EYE IMAGE CAPTU	JRING AND PROCESSI	NG DEVIC	E			
Application Infor	mation					
APPLICATION TYPE	Utility - Nonprovisional Application under 35 USC 111(a)		PATENT #			
CONFIRMATION #	8287		FILED BY	SU YU		
PATENT CENTER #	65106515		FILING DATE	12/28/2023		
CUSTOMER #	65358		FIRST NAMED	CHENG-CHUN CHANG		
CORRESPONDENCE ADDRESS	2		AUTHORIZED BY	JUSTIN KING		
Documents			ΤΟΤΑ	L DOCUMENTS: 2		
DOCUMENT		PAGES	DESCRIPTION	SIZE (KB)		
poa1.pdf		1	Power of Attorney	2601 KB		
poa2.pdf		1	Power of Attorney	4643 KB		
Digest						

DCUMENT	MESSAGE DIGEST(SHA-512)
a1.pdf	D464BFD67367C68880B1BCC48D4A693E8C8FED44C9ED1809
	76393526FFF7408CF6F92FE198A49B2EE5B21DFFA3C0C0D48
	2C5365995B7C0471C91B38CDC32AEA4

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